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Manufact  
TRB Chem

## INSTRUCTIONS FOR USE

VISMED® gel – Sodium Hyaluronate from fermentation C  
**Hydrogel for sustained lubrication of the eyes.** Str

**Composition:** Each 1 ml of solution contains: sodium hyaluronate, sodium chloride, disodium-hydrogenphosphate, sodium citrate and water for injections. The solution is hypotonic.

**Intended use:** For sustained lubrication of the eyes in case of ocular fatigue and other minor complaints of no pathological nature, caused e.g. by dust, smoke, dry heat, air conditioning, extended computer use or ophthalmological screening methods.

**Contra-Indications:** Individual hypersensitivity to any of the ingredients.

**Interactions:** Do not use VISMED® gel at the same time as other eye drops since it may modify their effects.

**Directions for use:** Twist off tab and position the tip of the container over the eye to be treated: upon slight pressure, VISMED® gel will easily flow out. Recommended, place one or two drops of VISMED® gel in the eye as needed. After blinking, the solution will disperse and form a thin film on the surface of the eye. VISMED® gel may also be used during contact lens wear.

**Precautions for use:** Do not touch the tip of the opening of the eye with the tip of the container. Do not use VISMED® gel if the solution contains preservatives, any solution should be discarded. Store between 2 °C and 25 °C! Do not use if the expiration date indicated on the package. If discomfort persists while using VISMED® gel, consult a physician. Keep out of reach of children.

**Characteristics and mode of action:** VISMED® gel contains hyaluronate, a natural polymer, which is also present in the structures of the human eye. The special physical characteristics of hyaluronate confer to VISMED® gel its "viscoelastic" and water-binding properties. VISMED® gel provides a stable coating on the surface of the eye, which is only gradually removed by blinking. Therefore VISMED® gel combines long lasting effect with maximum comfort. VISMED® gel is particularly compatible due to its unique composition. It is free of preservatives.

**Presentation:**  
20 monodose units of 0.45 ml; 60 monodose units of 0.45 ml  
Last revision date: August 2006

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TRB Chemedica AG, D - 85540 Haar/München