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D

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Sterile, free from preservatives.

#### Composition:

Each 1ml of solution contains: sodium hyaluronate 1.8 mg, sodium chloride, potassium chloride, disodium-hydrogenphosphate, sodium citrate, magnesium chloride, calcium chloride and water for injections. The solution is hypotonic.

#### Intended use:

Sensation of dryness, and other minor complaints of no pathological significance, as well as burning and ocular fatigue induced, for example, by dust, smoke, dry heat, air conditioning, extended computer screen use or contact lens wear.

#### Contra-Indications:

Individual hypersensitivity to any constituent of the product.

#### Interactions:

Do not use VISMED® at the same time with medicines or products applied to the eye since it may modify their effects.

#### Directions for use:

Twist off tab. If not otherwise recommended, place one or two drops of VISMED® into the conjunctival sac of the eye as often as needed. After blinking, the solution will disperse and form a transparent and long lasting coating on the surface of the eye. VISMED® may also be used while wearing contact lenses (rigid or soft).

#### Precaution for use:

Do not touch the tip of the container and do not touch the surface of the eye with the tip of the container. Do not use VISMED® if the container is damaged. As VISMED® does not contain preservatives, any solution not used immediately after opening should be discarded.

Store below 25°C! Do not use VISMED® after the expiry date indicated on the package. If discomfort persists while using VISMED® or reappears after discontinuation of administration, consult a physician.

Keep out of reach of children!

#### Characteristics and mode of action:

VISMED® contains sodium hyaluronate, a natural polymer, which is also present in the structures of the human eye. The particular physical characteristics of sodium hyaluronate confer to VISMED® its "viscoelastic" and water retaining properties. VISMED® provides a stable coating on the surface of the eye which is only slowly eliminated by blinking. Therefore VISMED® combines long lasting relief with maximum comfort. VISMED® is well tolerated due to its unique composition. It is free from preservatives.

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