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**(GB) How to clean my teeth properly**

Join the first animated, practical instruction according to the "New Swiss School of Oral hygiene" for a correct use of toothbrush, floss and interdental brushes. Designed according to the findings of iTOP. Visit [www.curaprox.com](http://www.curaprox.com) for more information.

- Hold your toothbrush clearly angulated against the gums. One half of the bristles covers the gums, the other half covers the surface of the teeth.
- The brush is slowly guided along the tooth-row from the front to the back by applying little pressure and fine circular movements.
- After the careful cleaning of the outside of the first half of the upper jaw you proceed to the other half.
- Proceed to cleaning the lingual surface of the upper teeth as for the buccal surfaces of the lower teeth.
- After the buccal cleansing, lingual cleansing is on.
- The compact head of the CS 820 ideally suits the vertical access to difficult to reach lingual anterior areas.
- While brushing the rear surfaces of the last molars the advantage of the compact head becomes particularly evident.
- The compact head allows perfect lingual cleansing even in case of narrow jaw anatomy.
- Last but not least. Brush all occlusal surfaces, particularly important for children.

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