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### GB AUS CA How to clean my teeth properly

Join the first animated, practical instruction for the correct use of toothbrush, floss and interdental brush. Designed according to the findings of iTOP.

Visit [www.curaprox.com](http://www.curaprox.com) for more information.

1. Hold your CURAPROX Single Brush at the right angle to your tooth, leading the outer edge of the bristle field carefully against the gum line.
2. Hold your Single Brush at a sloping angle after moving sideways.
3. Insert the outer edge of the brush very carefully against the gum line.
4. Slowly guide your Single Brush along the gum line with very fine jolting vertical movements.
5. At all times the bristle field must be kept in a vertical position against the tooth you are cleaning.
6. The cleaning action should start on the front-side on one of the sides of your tooth and finish on the opposite side of the same tooth.
7. Continue cleaning the tooth on the inside and clean the entire lingual surface.
8. Carefully feel your way along the gum line on the inside and clean your tooth according to the front-side.
9. The cleaning action inside will again start on one of the sides of your tooth and finish on the opposite side of the same tooth.

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